


Fluted Qualitative Cellulose



Increased efficiency in prep time

Easier to use and more convenient than flat filters when used in a filter funnel, this laboratory staple reduces your prep time with sharply creased prefolded pattern.

Applications:

- General laboratory separation & clarification
- Soil analysis methods
- Very high loading capacity
- Identical to Tremont grade CFP2

Retention

8µm

Flow speed

Medium-slow

Surface

Smooth

Brightness

Bright white

Ordering Information

I.W. Tremont catalog #:	Diameter (cm)	Packaging
CFP2F-125	12.5	100/pk
CFP2F-150	15.0	100/pk
CFP2F-185	18.5	100/pk
CFP2F-240	24.0	100/pk



Product is neatly packaged in 10 poly band groups of 10 ready to use filters.

Competitive Equivalents:

Description	Tremont catalog #	Company	Grades	Catalog #
125 mm	CFP2F-125	Ahlstrom	513	5130-1250
		GE/Whatman	2V	1202-125
150 mm	CFP2F-150	Ahlstrom	513	5130-1500
		GE/Whatman	2V	1202-150
185 mm	CFP2F-185	Ahlstrom	513	5130-1850
		GE/Whatman	2V	1202-185
240 mm	CFP2F-240	Ahlstrom	513	5130-2400
		GE/Whatman	2V	1202-240

Trademarks are owned by their respective owners. Use of competitive trademarks are for reference only – no ownership or rights are implied.

Cross reference information is a compilation of the listed manufacturers equivalency charts and sales data, not actual test data.

Although reasonable attempt has been made to ensure equivalency between competitive products - individual testing is suggested for non-standard method applications.