

Fluted Qualitative Cellulose

Increased efficiency in prep time

Easier to use and more convenient than flat filters when used in a filter funnel, this laboratory staple reduces your prep time with sharply creased prefolded pattern.


Applications:

General laboratory separation & clarification

Soil analysis methods

Very high loading capacity

Identical to Tremont grade CFP2



Retention

8µm

Flow speed

Medium-slow

Surface

Smooth

Brightness

Bright white

Ordering Information

I.W. Tremont catalog #:	Diameter (cm)	Packaging
CFP2F-125	12.5	100/pk
CFP2F-150	15.0	100/pk
CFP2F-185	18.5	100/pk
CFP2F-240	24.0	100/pk

Competitive Equivalents:

Description	Tremont catalog #	Company	Grades	Catalog #
125 mm	CFP2F-125	Ahlstrom	513	5130-1250
		GE/Whatman	2V	1202-125
150 mm	CFP2F-150	Ahlstrom	513	5130-1500
		GE/Whatman	2V	1202-150
185 mm	CFP2F-185	Ahlstrom	513	5130-1850
		GE/Whatman	2V	1202-185
240 mm	CFP2F-240	Ahlstrom	513	5130-2400
		GE/Whatman	2V	1202-240



Product is neatly packaged in 10 poly band groups of 10 ready to use filters.

Trademarks are owned by their respective owners. Use of competitive trademarks are for reference only – no ownership or rights are implied.

Cross reference information is a compilation of the listed manufacturers equivalency charts and sales data, not actual test data.

Although reasonable attempt has been made to ensure equivalency between competitive products - individual testing is suggested for non-standard method applications.